

# CDPF1 Antibody



PACO37138

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CDPF1

**Uniprot**

Q6NVV7

**Synonyms:**

Cysteine-rich DPF motif domain-containing protein 1, CDPF1, C22orf40

**Immunogen:**

Recombinant Human Cysteine-rich DPF motif domain-containing protein 1 protein (1-123AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A