

TRAT1 Antibody, HRP conjugated



PACO37163

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Stabilizes the TCR (T-cell antigen receptor)/CD3 complex at the surface of T-cells.

Gene ID:

TRAT1

Uniprot

Q6PIZ9

Synonyms:

T-cell receptor-associated transmembrane adapter 1 (T-cell receptor-interacting molecule) (TRIM) (pp29/30), TRAT1, TCRIM

Immunogen:

Recombinant Human T-cell receptor-associated transmembrane adapter 1 protein (29-186AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A