

SPATA8 Antibody



PACO37166

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

SPATA8

Uniprot

Q6RVD6

Synonyms:

Spermatogenesis-associated protein 8 (Spermatogenesis-related protein 8), SPATA8

Immunogen:

Recombinant Human Spermatogenesis-associated protein 8 protein (1-105AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A