

# C6orf89 Antibody



PACO37170

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Exhibits histone deacetylase (HDAC) enhancer properties. May play a role in cell cycle progression and wound repair of bronchial epithelial cells.

**Gene ID:**

C6orf89

**Uniprot**

Q6UWU4

**Synonyms:**

Bombesin receptor-activated protein C6orf89 (Amfion), C6orf89, BRAP

**Immunogen:**


Recombinant Human Bombesin receptor-activated protein C6orf89 protein (120-320AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human colon cancer using PACO37170 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO37170 at dilution of 1:100.