

DEF8 Antibody



PACO37178

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

DEF8

Uniprot

Q6ZN54

Synonyms:

Differentially expressed in FDCP 8 homolog (DEF-8), DEF8

Immunogen:

Recombinant Human Differentially expressed in FDCP 8 homolog protein (1-197AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A