

CEP68 Antibody, Biotin conjugated



PACO37201

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in maintenance of centrosome cohesion, probably as part of a linker structure which prevents centrosome splitting. Required for localization of CDK5RAP2 to the centrosome during interphase.

Gene ID:

CEP68

Uniprot

Q76N32

Synonyms:

Centrosomal protein of 68 kDa (Cep68), CEP68, KIAA0582

Immunogen:

Recombinant Human Centrosomal protein of 68 kDa protein (280-480AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A