

DPPA4 Antibody, Biotin conjugated



PACO37205

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in the maintenance of active epigenetic status of target genes. May inhibit differentiation of embryonic cells into a primitive ectoderm lineage.

Gene ID:

DPPA4

Uniprot

Q7L190

Synonyms:

Developmental pluripotency-associated protein 4, DPPA4

Immunogen:

Recombinant Human Developmental pluripotency-associated protein 4 protein (11-264AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A