

# BCDIN3D Antibody, Biotin conjugated



PACO37217

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

O-methyltransferase that specifically dimethylates the 5' monophosphate of pre-miRNAs, acting as a negative regulator of miRNA processing. The 5' monophosphate of pre-miRNAs is recognized by DICER1 and is required for pre-miRNAs processing; methylation at this position reduces the processing of pre-miRNAs by DICER1. Able to mediate methylation of pre-miR-145, as well as other pre-miRNAs.

**Gene ID:**

BCDIN3D

**Uniprot**

Q7Z5W3

**Synonyms:**

Pre-miRNA 5'-monophosphate methyltransferase (EC 2.1.1. -) (BCDIN3 domain-containing protein), BCDIN3D

**Immunogen:**

Recombinant Human Pre-miRNA 5'-monophosphate methyltransferase protein (1-292AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A