

KRT77 Antibody



PACO37218

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

KRT77

Uniprot

Q7Z794

Synonyms:

Keratin, type II cytoskeletal 1b (Cytokeratin-1B) (CK-1B) (Keratin-77) (K77) (Type-II keratin Kb39), KRT77, KRT1B


Immunogen:

Recombinant Human Keratin, type II cytoskeletal 1b protein (1-578AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human skin tissue using PACO37218 at dilution of 1:100.