

CPLX4 Antibody, HRP conjugated



PACO37223

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Positively regulates a late step in synaptic vesicle exocytosis.

Gene ID:

CPLX4

Uniprot

Q7Z7G2

Synonyms:

Complexin-4 (Complexin IV) (CPX IV), CPLX4

Immunogen:

Recombinant Human Complexin-4 protein (1-157AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A