

TADA2B Antibody, Biotin conjugated



PACO37233

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Coactivates PAX5-dependent transcription together with either SMARCA4 or GCN5L2.

Gene ID:

TADA2B

Uniprot

Q86TJ2

Synonyms:

Transcriptional adapter 2-beta (ADA2-like protein beta) (ADA2-beta), TADA2B, ADA2B

Immunogen:

Recombinant Human Transcriptional adapter 2-β protein (29-328AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A