

TSPAN33 Antibody, Biotin conjugated



PACO37237

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plays an important role in normal erythropoiesis. It has a role in the differentiation of erythroid progenitors.

Gene ID:

TSPAN33

Uniprot

Q86UF1

Synonyms:

Tetraspanin-33 (Tspan-33) (Penumbra) (hPen) (Proerythroblast new membrane), TSPAN33, PEN

Immunogen:

Recombinant Human Tetraspanin-33 protein (118-235AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A