

FAM122B Antibody



PACO37266

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:**Gene ID:**

FAM122B

Uniprot

Q7Z309

Synonyms:

Protein FAM122B, FAM122B


Immunogen:

Recombinant Human Protein FAM122B protein (2-247AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: FAM122B antibody at 5 μ g/ml + Mouse brain tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 27, 28, 30, 20 kDa. Observed band size: 27 kDa.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO37266 at dilution of 1:100.