

# ANKRD33 Antibody



PACO37274

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Acts as a transcriptional repressor for CRX-activated photoreceptor gene regulation.

**Gene ID:**

ANKRD33

**Uniprot**

Q7Z3H0

**Synonyms:**

Photoreceptor ankyrin repeat protein (Ankyrin repeat domain-containing protein 33), ANKRD33, C12orf7 PANKY

**Immunogen:**


Recombinant Human Photoreceptor ankyrin repeat protein (1-272AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO37274 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO37274 at dilution of 1:100.