

# CYP2U1 Antibody



PACO37278

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IF

**Recommended dilutions:**

ELISA:1:2000-1:10000, IF:1:50-1:200

**Protein Background:**

Catalyzes the hydroxylation of arachidonic acid, docosahexaenoic acid, and other long chain fatty acid. . May modulate the arachidonic acid, signaling pathway and play a role in other fatty acid, signaling processes.

**Gene ID:**

CYP2U1

**Uniprot**

Q7Z449

**Synonyms:**

Cytochrome P450 2U1 (EC 1.14.14.1), CYP2U1

**Immunogen:**


Recombinant Human Cytochrome P450 2U1 protein (363-496AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunofluorescent analysis of HeLa cells using PACO37278 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).