

TBC1D26 Antibody



PACO37298

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

May act as a GTPase-activating protein for Rab family protein(s).

Gene ID:

TBC1D26

Uniprot

Q86UD7

Synonyms:

TBC1 domain family member 26, TBC1D26


Immunogen:

Recombinant Human TBC1 domain family member 26 protein (1-250AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: TBC1D26 antibody at 3 μ g/ml + Mouse liver tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 29, 25, 10 kDa. Observed band size: 29 kDa.