

# VASH2 Antibody, HRP conjugated



PACO37303

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Angiogenesis inhibitor. Inhibits network formation by endothelial cells.

**Gene ID:**

VASH2

**Uniprot**

Q86V25

**Synonyms:**

Vasohibin-2 (Vasohibin-like protein), VASH2, VASHL

**Immunogen:**

Recombinant Human Vasohibin-2 protein (1-156AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A