

VASH2 Antibody, Biotin conjugated



PACO37305

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Angiogenesis inhibitor. Inhibits network formation by endothelial cells.

Gene ID:

VASH2

Uniprot

Q86V25

Synonyms:

Vasohibin-2 (Vasohibin-like protein), VASH2, VASHL

Immunogen:

Recombinant Human Vasohibin-2 protein (1-156AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A