

# SKAP1 Antibody, Biotin conjugated



PACO37317

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SKAP1

**Uniprot**

Q86WV1

**Synonyms:**

Src kinase-associated phosphoprotein 1 (Src family-associated phosphoprotein 1) (Src kinase-associated phosphoprotein of 55 kDa) (SKAP-55) (pp55), SKAP1, SCAP1 SKAP55

**Immunogen:**

Recombinant Human Src kinase-associated phosphoprotein 1 protein (1-358AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A