

# NUDT13 Antibody



PACO37318

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NUDT13

**Uniprot**

Q86X67

**Synonyms:**

Nucleoside diphosphate-linked moiety X motif 13 (Nudix motif 13) (EC 3. -. -. ) (Protein KiSS-16), NUDT13

**Immunogen:**

Recombinant Human Nucleoside diphosphate-linked moiety X motif 13 protein (32-263AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A