

# CSMD2 Antibody, HRP conjugated



PACO37327

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CSMD2

**Uniprot**

Q7Z408

**Synonyms:**

CUB and sushi domain-containing protein 2 (CUB and sushi multiple domains protein 2), CSMD2, KIAA1884

**Immunogen:**

Recombinant Human CUB and sushi domain-containing protein 2 protein (1-225AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A