

SPATA19 Antibody, HRP conjugated



PACO37339

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May have a role in spermiogenesis.

Gene ID:

SPATA19

Uniprot

Q7Z5L4

Synonyms:

Spermatogenesis-associated protein 19, mitochondrial (Spermatogenic cell-specific gene 1 protein) (Spergen-1), SPATA19, SPERGEN1

Immunogen:

Recombinant Human Spermatogenesis-associated protein 19, mitochondrial protein (25-167AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A