

# SPATA19 Antibody, FITC conjugated



PACO37340

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May have a role in spermiogenesis.

**Gene ID:**

SPATA19

**Uniprot**

Q7Z5L4

**Synonyms:**

Spermatogenesis-associated protein 19, mitochondrial (Spermatogenic cell-specific gene 1 protein) (Spergen-1), SPATA19, SPERGEN1

**Immunogen:**

Recombinant Human Spermatogenesis-associated protein 19, mitochondrial protein (25-167AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A