

PACO37346

Product Information

Size:

50ug

Reactivity:

Bacillus anthracis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plays an important role in DNA replication, recombination and repair. Binds to ssDNA and to an array of partner proteins to recruit them to their sites of action during DNA metabolism.

Gene ID:

ssb

Uniprot

Q81J13

Synonyms:

Single-stranded DNA-binding protein (SSB), ssb, ssb-1

Immunogen:

Recombinant Bacillus anthracis Single-stranded DNA-binding protein (1-172AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A