

PACO37348

---

## Product Information

**Size:**

50ug

**Reactivity:**

Bacillus anthracis

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Plays an important role in DNA replication, recombination and repair. Binds to ssDNA and to an array of partner proteins to recruit them to their sites of action during DNA metabolism.

**Gene ID:**

ssb

**Uniprot**

Q81J13

**Synonyms:**

Single-stranded DNA-binding protein (SSB), ssb, ssb-1

**Immunogen:**

Recombinant Bacillus anthracis Single-stranded DNA-binding protein (1-172AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A