

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Cytoplasmic protein involved in the delivery of the dynein machinery to the motile cilium. It is required for the assembly of the axonemal dynein inner and outer arms, two structures attached to the peripheral outer doublet A microtubule of the axoneme, that play a crucial role in cilium motility.

Gene ID:

DNAAF5

Uniprot

Q86Y56

Synonyms:

Dynein assembly factor 5, axonemal (HEAT repeat-containing protein 2), DNAAF5, HEATR2

Immunogen:

Recombinant Human Dynein assembly factor 5, axonemal protein (621-843AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A