

P4HA3 Antibody, HRP conjugated



PACO37387

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the post-translational formation of 4-hydroxyproline in -Xaa-Pro-Gly-sequences in collagens and other proteins.

Gene ID:

P4HA3

Uniprot

Q7Z4N8

Synonyms:

Prolyl 4-hydroxylase subunit alpha-3 (4-PH alpha-3) (EC 1.14.11.2) (Procollagen-proline,2-oxoglutarate-4-dioxygenase subunit alpha-3), P4HA3

Immunogen:

Recombinant Human Prolyl 4-hydroxylase subunit alpha-3 protein (20-544AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A