

GPR141 Antibody, HRP conjugated



PACO37451

Product Information

Size:	Protein Background:
50ug	Orphan receptor.
Reactivity:	Gene ID:
Human	GPR141
Source:	Uniprot
Rabbit	Q7Z602
Isotype:	Synonyms:
IgG	Probable G-protein coupled receptor 141 (G-protein coupled receptor PGR13), GPR141, PGR13
Applications:	Immunogen:
ELISA	Recombinant Human Probable G-protein coupled receptor 141 protein (153-183AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A