

# GPR141 Antibody, Biotin conjugated



PACO37453

---

## Product Information

<b>Size:</b>	<b>Protein Background:</b>
50ug	Orphan receptor.
<b>Reactivity:</b>	<b>Gene ID:</b>
Human	GPR141
<b>Source:</b>	<b>Uniprot</b>
Rabbit	Q7Z602
<b>Isotype:</b>	<b>Synonyms:</b>
IgG	Probable G-protein coupled receptor 141 (G-protein coupled receptor PGR13), GPR141, PGR13
<b>Applications:</b>	<b>Immunogen:</b>
ELISA	Recombinant Human Probable G-protein coupled receptor 141 protein (153-183AA).
<b>Recommended dilutions:</b>	<b>Storage:</b>
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A