

LHFPL1 Antibody



PACO37470

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Gene ID:

LHFPL1

Uniprot

Q86WI0

Synonyms:

LHFPL tetraspan subfamily member 1 protein (Lipoma HMGIC fusion partner-like 1 protein), LHFPL1


Immunogen:

Recombinant Human LHFPL tetraspan subfamily member 1 protein (21-85AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Immunofluorescent analysis of HeLa cells using PACO37470 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human liver cancer using PACO37470 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO37470 at dilution of 1:100.