

PID1 Antibody, FITC conjugated



PACO37492

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Increases proliferation of preadipocytes without affecting adipocytic differentiation.

Gene ID:

PID1

Uniprot

Q7Z2X4

Synonyms:

PTB-containing, cubilin and LRP1-interacting protein (P-CLI1) (Phosphotyrosine interaction domain-containing protein 1) (Protein NYGGF4), PID1, NYGGF4 PCLI1

Immunogen:

Recombinant Human PTB-containing, cubilin and LRP1-interacting protein (15-168AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A