

RD3 Antibody, Biotin conjugated



PACO37501

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

RD3

Uniprot

Q7Z3Z2

Synonyms:

Protein RD3 (Retinal degeneration protein 3), RD3, C1orf36

Immunogen:

Recombinant Human Protein RD3 protein (1-195AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A