

VMO1 Antibody, FITC conjugated



PACO37504

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

VMO1

Uniprot

Q7Z5L0

Synonyms:

Vitelline membrane outer layer protein 1 homolog, VMO1

Immunogen:

Recombinant Human Vitelline membrane outer layer protein 1 homolog protein (25-202AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A