

# YTHDF3 Antibody, HRP conjugated



PACO37507

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Specifically recognizes and binds N6-methyladenosine (m6A)-containing RNAs. M6A is a modification present at internal sites of mRNAs and some non-coding RNAs and plays a role in the efficiency of mRNA splicing, processing and stability.

**Gene ID:**

YTHDF3

**Uniprot**

Q7Z739

**Synonyms:**

YTH domain-containing family protein 3, YTHDF3

**Immunogen:**

Recombinant Human YTH domain-containing family protein 3 protein (315-585AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A