

YTHDF3 Antibody, FITC conjugated



PACO37508

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Specifically recognizes and binds N6-methyladenosine (m6A)-containing RNAs. M6A is a modification present at internal sites of mRNAs and some non-coding RNAs and plays a role in the efficiency of mRNA splicing, processing and stability.

Gene ID:

YTHDF3

Uniprot

Q7Z739

Synonyms:

YTH domain-containing family protein 3, YTHDF3

Immunogen:

Recombinant Human YTH domain-containing family protein 3 protein (315-585AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A