

TRIM74 Antibody



PACO37514

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

TRIM74

Uniprot

Q86UV6

Synonyms:

Tripartite motif-containing protein 74 (Tripartite motif-containing protein 50C), TRIM74, TRIM50C


Immunogen:

Recombinant Human Tripartite motif-containing protein 74 protein (47-249AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO37514 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO37514 at dilution of 1:100.