

MON1A Antibody, HRP conjugated



PACO37523

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays an important in membrane trafficking through the secretory apparatus. Not involved in endocytic trafficking to lysosomes.

Gene ID:

MON1A

Uniprot

Q86VX9

Synonyms:

Vacuolar fusion protein MON1 homolog A, MON1A, SAND1

Immunogen:

Recombinant Human Vacuolar fusion protein MON1 homolog A protein (205-555AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A