

# MON1A Antibody, Biotin conjugated



PACO37525

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Plays an important in membrane trafficking through the secretory apparatus. Not involved in endocytic trafficking to lysosomes.

**Gene ID:**

MON1A

**Uniprot**

Q86VX9

**Synonyms:**

Vacuolar fusion protein MON1 homolog A, MON1A, SAND1

**Immunogen:**

Recombinant Human Vacuolar fusion protein MON1 homolog A protein (205-555AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A