

CACNA2D4 Antibody, Biotin conjugated



PACO37545

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

The alpha-2/delta subunit of voltage-dependent calcium channels regulates calcium current density and activation/inactivation kinetics of the calcium channel.

Gene ID:

CACNA2D4

Uniprot

Q7Z3S7

Synonyms:

Voltage-dependent calcium channel subunit alpha-2/delta-4 (Voltage-gated calcium channel subunit alpha-2/delta-4) [Cleaved into: Voltage-dependent calcium channel subunit alpha-2-4; Voltage-dependent calcium channel subunit delta-4], CACNA2D4

Immunogen:

Recombinant Human Voltage-dependent calcium channel subunit alpha-2/delta-4 protein (934-1115AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A