

PACO37550

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Involved in the regulation of membrane traffic through the trans-Golgi network (TGN).

Gene ID:

CCDC91

Uniprot

Q7Z6B0

Synonyms:

Coiled-coil domain-containing protein 91 (GGA-binding partner) (p56 accessory protein), CCDC91, GGABP


Immunogen:

Recombinant Human Coiled-coil domain-containing protein 91 protein (1-180AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human ovarian cancer using PACO37550 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human melanoma using PACO37550 at dilution of 1:100.