

PACO37566

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:1000-1:5000,  
IHC:1:20-1:200

**Protein Background:**

Required for polyglutamylation of axonemal tubulin. Plays a role in anterograde intraflagellar transport (IFT), the process by which cilia precursors are transported from the base of the cilium to the site of their incorporation at the tip.

**Gene ID:**

TTC30A

**Uniprot**

Q86WT1

**Synonyms:**

Tetratricopeptide repeat protein 30A (TPR repeat protein 30A), TTC30A

**Immunogen:**


Recombinant Human Tetratricopeptide repeat protein 30A protein (201-500AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images


---



Western blot. All lanes: TTC30A antibody at 0.1 $\mu$ g/ml + Mouse gonadal tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 77 kDa. Observed band size: 77 kDa.



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO37566 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO37566 at dilution of 1:100.