

Paox Antibody, HRP conjugated



PACO37583

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Flavoenzyme which catalyzes the oxidation of N(1)-acetylspermine to spermidine and is thus involved in the polyamine back-conversion. Can also oxidize N(1)-acetylspermidine to putrescine. Substrate specificity: N(1)-acetylspermine = N(1)-acetylspermidine > N(1), N(12)-diacylspermine >> spermine. Does not oxidize spermidine. Plays an important role in the regulation of polyamine intracellular concentration and has the potential to act as a determinant of cellular sensitivity to the antitumor polyamine analogs.

Gene ID:

Paox

Uniprot

Q8C0L6

Synonyms:

Peroxisomal N(1)-acetyl-spermine/spermidine oxidase (EC 1.5.3.13) (Polyamine oxidase), Paox, Pao

Immunogen:

Recombinant Mouse Peroxisomal N(1)-acetyl-spermine/spermidine oxidase protein (1-504AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A