

TMPRSS6 Antibody, HRP conjugated



PACO37591

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Serine protease which hydrolyzes a range of proteins including type I collagen, fibronectin and fibrinogen. Can also activate urokinase-type plasminogen activator with low efficiency. May play a specialized role in matrix remodeling processes in liver. Through the cleavage of HFE2, a regulator of the expression of the iron absorption-regulating hormone hepcidin/HAMP, plays a role in iron homeostasis.

Gene ID:

TMPRSS6

Uniprot

Q8IU80

Synonyms:

Transmembrane protease serine 6 (EC 3.4.21. -) (Matriptase-2), TMPRSS6

Immunogen:

Recombinant Human Transmembrane protease serine 6 protein (98-301AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A