

C1orf210 Antibody



PACO37602

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in membrane trafficking between endosomes and plasma membrane.

Gene ID:

C1orf210

Uniprot

Q8IVY1

Synonyms:

Type III endosome membrane protein TEMP (TEMP), C1orf210

Immunogen:

Recombinant Human Type III endosome membrane protein TEMP protein (51-113AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A