

ADGRA3 Antibody, FITC conjugated



PACO37608

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Orphan receptor that may have a role in planar cell polarity pathway.

Gene ID:

ADGRA3

Uniprot

Q8IWK6

Synonyms:

Adhesion G protein-coupled receptor A3 (G-protein coupled receptor 125), ADGRA3, GPR125

Immunogen:

Recombinant Human Adhesion G protein-coupled receptor A3 protein (101-381AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A