

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Key component of the nexin-dynein regulatory complex (N-DRC), essential for N-DRC integrity. Involved in the regulation of flagellar motility.

Gene ID:

DRC7

Uniprot

Q8IY82

Synonyms:

Dynein regulatory complex subunit 7 (Coiled-coil domain-containing protein 135) (Coiled-coil domain-containing protein lobo homolog), DRC7, C16orf50 CCDC135


Immunogen:

Recombinant Human Dynein regulatory complex subunit 7 protein (1-300AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO37626 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO37626 at dilution of 1:100.