

PACO37682

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

May play an important role in lymphopoiesis.

**Gene ID:**

RTKN2

**Uniprot**

Q8IZC4

**Synonyms:**

Rhotekin-2 (Pleckstrin homology domain-containing family K member 1) (PH domain-containing family K member 1), RTKN2, PLEKHK1

**Immunogen:**


Recombinant Human Rhotekin-2 protein (1-163AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human testis tissue using PACO37682 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human skin tissue using PACO37682 at dilution of 1:100.