

# EGFL6 Antibody, HRP conjugated



PACO37691

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May bind integrin alpha-8/beta-1 and play a role in hair follicle morphogenesis. Promotes matrix assembly.

**Gene ID:**

EGFL6

**Uniprot**

Q8IUJ8

**Synonyms:**

Epidermal growth factor-like protein 6 (EGF-like protein 6) (MAM and EGF domains-containing gene protein), EGFL6, MAEG

**Immunogen:**

Recombinant Human Epidermal growth factor-like protein 6 protein (301-554AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A