

EGFL6 Antibody, FITC conjugated



PACO37692

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May bind integrin alpha-8/beta-1 and play a role in hair follicle morphogenesis. Promotes matrix assembly.

Gene ID:

EGFL6

Uniprot

Q8IUJ8

Synonyms:

Epidermal growth factor-like protein 6 (EGF-like protein 6) (MAM and EGF domains-containing gene protein), EGFL6, MAEG

Immunogen:

Recombinant Human Epidermal growth factor-like protein 6 protein (301-554AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A