

NUDCD3 Antibody



PACO37694

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

NUDCD3

Uniprot

Q8IVD9

Synonyms:

NudC domain-containing protein 3, NUDCD3, KIAA1068


Immunogen:

Recombinant Human NudC domain-containing protein 3 protein (142-361AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: NUDCD3 antibody at 5 μ g/ml. Lane 1: 293T whole cell lysate. Lane 2: Hela whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 41 kDa. Observed band size: 41 kDa.